


Disinfectant Spray (for health care use) S3324

A phenol-based, pleasantly scented product for disinfecting after cleaning. It effectively disinfects HIV-1 (AIDS Virus), Pandemic 2009 H1N1 Influenza A, Influenza A2/Hong Kong Virus and Herpes Simplex Virus Types 1 and 2, Infectious Canine Hepatitis, Canine Distemper, Avian Influenza A (H5N1) on inanimate environmental surfaces. It is also Tuberculocidal, bactericidal, staphylocidal, fungicidal, pseudomonicidal, air deodorizer and is effective in controlling mold and mildew on hard, nonporous surfaces.



Where To Use

This product disinfects and deodorizes as you spray. This EPA registered product is great for use in/on:

Schools	Telephones	Locker Rooms	Athletic Equipment
Kitchens	Bathrooms	Garbage Cans	Day Care Centers
Hospitals	Handles	Toilet Seats	Shower Curtains
Offices			

Directions

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that is introduced directly into the human body.

Product is ready to use. Shake container well before each use and hold it upright while spraying. For use on precleaned, hard nonporous surfaces; not for use on humans or animals. Preclean surfaces Hold can about 6-8" from surface. Spray to thoroughly wet hard, nonporous surfaces. Treated surfaces must remain wet for 10 minutes. Allow to air dry. Rinse all surfaces that come in contact with food such as counters, appliances, tables with potable water. Do not use on utensils, glassware, and dishes.

Technical Specifications

CAN SIZE: ...20 oz.	NET WEIGHT: 16.5 oz.
ODOR:Pleasant	USDA RATING: C1
TYPE:Aerosol	PROPELLANT:Hydrocarbon
pH:7.0	SPRAY PATTERN: .Misty Solid Cone
COLOR:Clear	FLAMMABILITY:.... Flammable - as described in 40 CFR 156.10



201 Boston Turnpike, Bolton, CT 06043-7203 800-227-5536 • 860-646-0172 • FAX 860-645-6070

www.simoniz.com